

CHNS-ANALYSIS / O-ANALYSIS

References:

Instruction manual EA 1110

Elemental Analyzers (CE Instruments).

CHNS-Analysis

The determination of Carbon, Hydrogen, Nitrogen and Sulfur content is carried out by a simultaneous analysis.

- 1 Solid and liquid samples, weighed in tin capsules, are introduced into a vertical quartz reactor heated at a temperature of 1020°C with a constant flow of helium stream.
- 2 A few seconds before introduction the helium stream was enriched with high purity oxygen.
- 3 The combustion gas mixture is driven through an tungsten oxide zone to achieved a complete quantitative oxidation following by a reduction step in a copper

zone to reduce nitrogen oxides and sulfuric anhydride to nitrogen and sulfurous anhydride.

- 4 The resulting four components N_2 , CO_2 , H_2O und SO_2 are separated in a chromatographic column and detected by a thermo conductivity detector.
- 5 The resulting signals, proportional to the amount of eluted gases, are analyzed by an automatic workstation which provides the sample elemental composition report.

O-Analysis

The determination of Oxygen is a separate analysis.

- 1 The samples are weighed in silver capsules.
- 2 These capsules are dropped into a pyrolysis reactor kept at a temperature of 1060°C and crossed by a helium stream.
- 3 The samples are undergone an immediate pyrolysis
- 4 Pyrolysis gases are passed on a nickel-plated carbon layer that ensures a quantitative conversion of organic oxygen into carbon monoxide.
- 5 Carbon monoxide is separated from the other pyrolysis gases (CH_4 , N_2 , H_2 and other acid gases) in a GC column.
- 6 The gases are detected by a thermo conductivity detector.
- 7 The resulting signals, proportional to the amount of eluted gases, are analyzed by an automatic workstation which provides the sample elemental composition report.

Fraunhofer Institute for Applied Polymer Research IAP

Science Park Potsdam-Golm
Geiselbergstr. 69
14476 Potsdam-Golm

Contact

Dr. Hendrik Wetzel

Phone +49 331 568-1604
hendrik.wetzel@iap.fraunhofer.de

www.iap.fraunhofer.com




Fig. 1 CHNS-Analysis