



## MOLAR MASS DETERMINATION FROM SEDIMENTATION DIFFUSION EQUILIBRIUM

With moderate rotor speed at sufficiently long run time of the centrifuge an equilibrium between sedimentation and back-diffusion of the dissolved macromolecules is reached. One distinguishes between "high speed" or "meniscus depletion" equilibrium, where the polymer concentration at the meniscus comes to zero, and "low speed" equilibrium (Fig. 1) which is more adequate for samples with broad mass distribution. From the equilibrium concentration distribution between cell meniscus and bottom which is detected by UV-absorbance or interference optics the average molar masses  $M_w$  (weight average) and  $M_z$  (centrifuge average) can be calculated after an appropriate extrapolation to zero concentration (Fig. 2).



### Fraunhofer Institute for Applied Polymer Research IAP

Science Park Potsdam-Golm  
Geiselbergstr. 69  
14476 Potsdam-Golm

Contact

#### Dr. Erik Wischerhoff

Phone +49 331 568-1508  
erik.wischerhoff@iap.fraunhofer.de

[www.iap.fraunhofer.com](http://www.iap.fraunhofer.com)